

## EnzyFluo™ NADP<sup>+</sup>/NADPH Assay Kit (EFNP-100)

### Quantitative Fluorimetric Determination of NADP<sup>+</sup>/NADPH

#### DESCRIPTION

Pyridine nucleotides play an important role in metabolism and, thus, there is continual interest in monitoring their concentration levels. Quantitative determination of NADP<sup>+</sup>/NADPH has applications in research pertaining to energy transformation and redox state of cells or tissue.

Simple, direct and automation-ready procedures for measuring NADP<sup>+</sup>/NADPH concentration are very desirable. BioAssay Systems' EnzyFluo™ NADP<sup>+</sup>/NADPH assay kit is based on a glucose dehydrogenase cycling reaction, in which the formed NADPH reduces a probe into a highly fluorescent product. The fluorescence intensity of this product, measured at  $\lambda_{\text{ex/em}} = 530/585$  nm, is proportional to the NADP<sup>+</sup>/NADPH concentration in the sample. This assay is highly specific for NADP<sup>+</sup>/NADPH and with minimal interference (<1%) by NAD<sup>+</sup>/NADH. Our assay is a convenient method to measure NADP<sup>+</sup>, NADPH and their ratio.

#### APPLICATIONS

**Direct Assays:** NADP<sup>+</sup>/NADPH concentrations and ratios in cell or tissue extracts.

#### KEY FEATURES

**Sensitive and accurate.** Detection limit of 0.01  $\mu\text{M}$  and linearity up to 1  $\mu\text{M}$  NADP<sup>+</sup>/NADPH in 96-well plate assay.

**Convenient.** The procedure involves adding a single working reagent, and reading the fluorescence at time zero and 30 min. Room temperature assay.

**High-throughput.** Can be readily automated as a high-throughput 96-well plate assay for thousands of samples per day.

#### KIT CONTENTS

**Assay Buffer:** 10 mL      **Enzyme A:** 120  $\mu\text{L}$   
**Glucose (1 M):** 1.5 mL      **Enzyme B:** 120  $\mu\text{L}$   
**Probe:** 750  $\mu\text{L}$       **NADP Standard:** 0.5 mL  
**NADP/NADPH Extraction Buffers:** each 12 mL

**Storage conditions.** The kit is shipped on ice. Store all reagents at -20°C. Shelf life: 6 months after receipt.

**Precautions:** reagents are for research use only. Normal precautions for laboratory reagents should be exercised while using the reagents. Please refer to Material Safety Data Sheet for detailed information.

#### GENERAL CONSIDERATIONS

- At these concentrations, the standard curves for NADP<sup>+</sup> and NADPH are identical. Since NADPH in solution is unstable, we provide only NADP as the standard.
- This assay is based on an enzyme-catalyzed kinetic reaction. Addition of Working Reagent should be quick and mixing should be brief but thorough. Use of multi-channel pipettor is recommended.
- The following substances interfere and should be avoided in sample preparation. EDTA (>0.5 mM), ascorbic acid, SDS (>0.2%), sodium azide, NP-40 (>1%) and Tween-20 (>1%).

#### PROCEDURES

- Sample Preparation.** For tissues weigh ~20 mg tissue for each sample, wash with cold PBS. For cell samples, wash cells with cold PBS and pellet ~10<sup>5</sup> cells for each sample. Homogenize samples (either tissue or cells) in a 1.5 mL Eppendorf tube with either 100  $\mu\text{L}$  NADP<sup>+</sup> extraction buffer for NADP<sup>+</sup> determination or 100  $\mu\text{L}$  NADPH extraction buffer for NADPH determination. Heat extracts at 60°C for 5 min and then add 20  $\mu\text{L}$  Assay Buffer and 100  $\mu\text{L}$  of the opposite extraction buffer to neutralize the extracts. Briefly vortex and spin the samples down at 14,000 rpm for 5 min. Use supernatant for NADP<sup>+</sup>/NADPH assays. Determination of both NADP<sup>+</sup> and NADPH concentrations requires extractions from two separate samples.
- Calibration Curve.** Prepare 5000  $\mu\text{L}$  1  $\mu\text{M}$  NADP<sup>+</sup> Premix by mixing 5  $\mu\text{L}$  1 mM Standard and 4995  $\mu\text{L}$  distilled water. Dilute standard as follows.

No	Premix + H <sub>2</sub> O	NADP ( $\mu\text{M}$ )
1	100 $\mu\text{L}$ + 0 $\mu\text{L}$	1.0
2	60 $\mu\text{L}$ + 40 $\mu\text{L}$	0.6
3	30 $\mu\text{L}$ + 70 $\mu\text{L}$	0.3
4	0 $\mu\text{L}$ + 100 $\mu\text{L}$	0

Transfer 50  $\mu\text{L}$  standards into wells of a black flat-bottom 96-well plate.

- Samples.** Add 50  $\mu\text{L}$  sample per well in separate wells.
- Reagent Preparation.** For each reaction well, prepare Working Reagent by mixing 40  $\mu\text{L}$  Assay Buffer, 1  $\mu\text{L}$  Enzyme A, 1  $\mu\text{L}$  Enzyme B, 10  $\mu\text{L}$  Glucose and 5  $\mu\text{L}$  Probe. Fresh reconstitution is recommended.
- Reaction.** Add 50  $\mu\text{L}$  Working Reagent per well quickly. Tap plate to mix.
- Read fluorescence at  $\lambda_{\text{ex/em}} = 530/585$  nm for time "zero" ( $F_0$ ) and  $F_{30}$  after a 30-min incubation at room temperature. Protect plate from light during this incubation.

#### CALCULATION

First compute the  $\Delta F$  for each Subtract standard and sample by subtracting  $F_0$  from  $F_{30}$ . Plot the standard  $\Delta F$ 's and determine the slope. The NADP<sup>+</sup> (H) concentration of the sample is computed as follows:


$$[\text{NADP(H)}] = \frac{\Delta F_{\text{SAMPLE}} - \Delta F_{\text{BLANK}}}{\text{Slope } (\mu\text{M}^{-1})} \times n \quad (\mu\text{M})$$

where  $\Delta F_{\text{SAMPLE}}$  and  $\Delta F_{\text{BLANK}}$  are the change in fluorescence intensity values of the Sample and Blank (STD 4) respectively. Slope is the slope of the standard curve and  $n$  is the dilution factor (if necessary).

Note: If the sample  $\Delta F$  values are higher than the  $\Delta F$  value for the 1  $\mu\text{M}$  standard, dilute sample in distilled water and repeat this assay. Multiply the results by the dilution factor.

#### MATERIALS REQUIRED, BUT NOT PROVIDED

Pipetting (multi-channel) devices. Black, flat bottom 96-well plates and fluorescent plate reader capable of reading at  $\lambda_{\text{ex/em}} = 530/585$  nm.



Standard Curve in 96-well plate assay

#### PUBLICATIONS

- Li, Xiao-xue, et al (2018). Nuclear receptor Nur77 facilitates melanoma cell survival under metabolic stress by protecting fatty acid oxidation. *Molecular cell* 69.3: 480-492.
- Scherzi, T et al. (2021). Staphylococcus aureus resistance to albocycline can be achieved by mutations that alter cellular NAD/PH pools. *Bioorganic & Medicinal Chemistry*, 32, 115995.
- Guo, Shuang, et al. (2015) Resveratrol attenuates high glucose-induced oxidative stress and cardiomyocyte apoptosis through AMPK. *Molecular and cellular endocrinology* 412: 85-94.